



Thru-the-Lens Headlight Video Imaging System

The technological breakthrough
in surgical headlight imaging

Easy to use, our exclusive
Parallax Free System provides
an image of unsurpassed quality
in a single-chip headlight
video camera



Intelligently bright. Very bright.

BFW, INC. / 2307 RIVER ROAD, STE. 103 / LOUISVILLE, KY USA / 40206 WWW.BFWINC.COM / 502.899.1808 / 800.717.4673



Thru-the-Lens Headlight Video Imaging System

Available with tower cart for a flat screen monitor and recorder.

Our camera is custom engineered for use with a surgical headlight – not a refitted endoscopic or pre-built surveillance camera. This is how we offer incomparable single-chip image quality. With DV (IEEE 1394), S-VHS and BNC capability, our system offers an extended depth of focus. The image remains in focus regardless of inward or outward movement from the field of surgery.

The Thru-the-Lens Headlight Video Camera System is offered exclusively with the ChromaLUME™ TURBO Light Source with Plasma Lamp Technology and 10,000 hour lamp life.

S p e c i f i c a t i o n s

The Thru-the-Lens Headlight Video Camera system includes:

- Thru-the-Lens Headlight with 10' Fiber Optic Cable
- Camera Control Unit
- ChromaLUME™ TURBO Light Source
- All necessary Video Cables and Power Cords
- Mobile Stand

Product Codes

BFW 7930

Thru-the-Lens Headlight Video Camera System (NTSC)
With ChromaLUME™ TURBO Light Source

BFW 7935

Thru-the-Lens Headlight Video Camera System (PAL)
With ChromaLUME™ TURBO Light Source

MID 8610

10' Replacement Fiber Optic Cable

BFW 1505

Medical Grade DVD Recorder

BFW 1516

Flat Screen Monitor

BFW 1520

Tower Cart



Intelligently bright. Very bright.