



# VistaView II

## Portable LED Headlight System

- Lightweight & Portable
- Very Bright LED Illumination
- Sharp, Clean, Even Spot
- 2.2-watt LED
- 6000° Kelvin
- 50,000 hour lamp life
- 2 Rechargeable Batteries Included



*Intelligently bright. Very bright.*



# VistaView II

## Portable LED Headlight System

### Four Brightness Settings/Operating Times

- 18 hours - 9,000 Lux/900 foot candles
- 9 hours - 16,000 Lux/1,600 foot candles
- 7 hours - 20,000 Lux/2,000 foot candles
- 5 hours - 24,000 Lux/2,400 foot candles

*(ALL MEASUREMENTS TAKEN @ 16" (40 CM) WORKING DISTANCE)*

- Bright, White Light
- Sharp, Clean, Even Illumination
- Weight Battery: 125g, (4.4 oz)  
LED Optic: 35.1g (1.2 oz)
- Charging Time: 3.5 - 5 Hours
- Battery Life Span 300 Cycles
- Cable Length 1.5m/5ft
- 110/240 VAC, 50.60Hz
- Can Charge Battery While Light is in Use
- In-line Power and Intensity Control  
Within Easy Reach
- Warranty: Battery Case 1-year;  
Optic 2-years

### **BFW 9820**

System includes: LED optic, adjustable headband with brethe-o-prene™ comfort cushions, adapter to clip onto loupe glasses, 2 rechargeable battery packs, Battery Control Unit and in-line power and intensity control, charging stand, power supply for charging battery.

### **BFW 9825**

System includes: LED Optic, breathe-o-prene™ headband, adapter to clip onto loupe glasses, 2 rechargeable battery packs, Battery Control Unit and in-line power and intensity control, charging stand, power supply for charging battery.

### **BFW 9821**

As described above without headband.

### **BFW 9815**

Replacement Headband with comfort cushions.



*Intelligently bright. Very bright.*